


Live Webinar: 3D Tumor Spheroid based Functional Assays using Celigo Image Cytometer

Wednesday, October 21 | 2:00pm-3:00pm

Fast, Whole Well Imaging for 3D Tumor Spheroid Screening.




Celigo Imaging Cell Cytometer

Monitor spheroid growth over time.
Screen multiple spheroid plates per day.

Minimally Invasive, Microplate-Based Cytometry

- High resolution — 1 micron/pixel
- Every cell in every well — 96-well plate in < 5 minutes.
- 3 fluorescent colors plus bright field.
- Image 3D tumor spheroids, adherent and suspension cells.



Contate-nos para maiores informações!