


Sample System Overview MP160 with Wired Biopotential Amp



New features:

- ✓ Improved signal quality. Lower noise also improves signal quality for active electrodes.
- ✓ Includes **AcqKnowledge 5** software
- ✓ Auto-sensing Ethernet port. The MP160 can detect whether or not a patch cable or a crossover cable is connected automatically and work with either cable.
- ✓ **Includes an HLT100C** High Level Transducer Module
- ✓ Adjustable triggering hysteresis. This allows for finer adjustment of the window used for triggering detection off of analog input channels, similar to control calculation channels. In the MP150 this window is not adjustable.
- ✓ Lower power consumption.
- ✓ RoHS compliant/eco-friendly. Contains no lead.
- ✓ Expanded internal logic to allow new functionality to be added in the future.
- ✓ Communicates via Ethernet (USB adapter)

Contate-nos e peça já sua cotação!